

iTero Packaging Instructions



Required parts

Foam Base



Foam Shelf



Wood Pallet



Wand Box



Wand Thermal Cover



Scanner Foam Support



Tube Set



Monitor Package




External Carton Box & Cover





Insert foam base onto wooden surface.



Insert system into foam base.



Insert foam shelf support onto keyboard.





Insert scanner into the blanket's pocket



Tightly roll the wand in the blanket



Insert wand into box



Place wand on foam shelf according to the picture.






Place foot pedal in base at the front of the cart.



Insert the monitor into a plastic bag and place the original foam supports on monitor as shown.



Place monitor box over the screen.



Insert support poles





Insert carton around system as shown.



Assemble top cover.



Strap scanner to wood pallet

